MultiMotion

Dynamic Corrective Joints

Professional instructions

Basko
Healthcare


basko.com

MultiMotion

Dynamic Corrective Joints

Fig. A

1.



2.


3.

4.

5.

6.

NL

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DE

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ENG

Technical manual.....25

MultiMotion

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MultiMotion

Please read these instructions carefully before using the product. Pay particular attention that the safety instructions are followed.

Usage guidelines

- The product may only be fitted by or under the supervision of a competent professional
- Special precautions must always be taken for patients with a visual impairment, cognitive limitations and/or reduced sensibility in the upper and/or lower extremity.
- Changes, modifications and adjustments to the product that are not described in these instructions, are done at the sole responsibility of the person carrying them out.

Intended use of the product

By using MultiMotion dynamic correctable system joints – hereafter: corrective joints – dynamically correctable contractures in the limbs can be treated. The corrective joints are to be used exclusively in orthotics for the upper or lower extremities and should not be weight bearing. For correct functioning, it should always be combined with a MultiMotion free motion joint / R.O.M. joint at the contralateral side of the affected joint and ensures sufficient torsional stiffness of the orthosis.

Indications / contraindications and function

Multimotion joints are suitable for use in orthotics for wrist, elbow, hip, knee and ankle joints for adults and children. The joints are used for the treatment of dynamic correctable contractures resulting from neurological and/or orthopedic disorders.

Indications:

- | | |
|--------------------------------|----------------------|
| ■ Apoplexy (CVA) | (ICD10: I64) |
| ■ Burns | (ICD10: T95:9) |
| ■ Cerebral palsy | (ICD10: G80.9) |
| ■ Cranial Brain Trauma | (ICD10: S06.9) |
| ■ Endoprosthesis implant | |
| ■ Fracture | (ICD10: T14.2) |
| ■ Hereditäre spastische Parese | (ICD10: G11.4) |
| ■ Hip subluxation, congenital | (ICD10: Q65.5) |
| ■ Ligament rupture | (ICD10: T14.-) |
| ■ Multiple sclerosis | (ICD10: G35) |
| ■ Spastic diplegia | (ICD10: G80.1) |
| ■ Spastische Parese | (ICD10: G82.1) |
| ■ Spastische Tetraparese | (ICD10: G82.4) |
| ■ Spina Bifida | (ICD10: Q05) |
| ■ Spinal paralysis, spastic | (ICD10: G11.4) |
| ■ Tetraplegie, spastisch | (ICD10: G80.0) |
| ■ Transtibial amputation | (ICD10: Z89.5/Z89.7) |
| ■ Triplegie, spastisch | (ICD10: G80.9) |

The diagnosis is made by the physician or a competent medical professional.

Contraindications

- | | |
|---------------------------------|-----------------|
| ■ Ankylosis | (ICD10: M24.69) |
| ■ Fibrosis | (ICD10: M35) |
| ■ Muscles and bones deformation | |
| ■ Ossification | (ICD10: M61.9) |

General safety instructions

Risk of injury as a result of incorrect use; the professional should verify correct functioning of the joints prior to delivery of the orthosis to the user. The patient and his/her care giver(s) should also be instructed on donning and doffing of the orthosis and have the joint functionality explained to them. It should also be made clear to them that the spring tension setting should never be changed by anybody else than the professional or therapist, unless the treating physician decides otherwise. The joints should not come into contact with water.

Comment

Specialists opting for and/or fitting and customize these products, should use their professional judgment in choosing the product, fitting it and adequately educating the patient or care giver so as to minimize potential risks for each individual patient. These risks may include the mentioned contraindications, but may also be associated with the unique patient characteristics or circumstances of the patient's care giver.

Used materials

The MultiMotion functional corrective joints incorporate components made of various types of materials, including:

- Aluminum
- Steel/stainless steel
- Brass
- Synthetic PC ABS

Maintenance recommendation

The tension spring in the orthosis has already been greased during joint assembly. If the tension spring needs to be greased again, a standard white grease (e.g. Molykote) should be used. Insufficient grease on the tension spring may cause undesired friction sounds.


We recommend checking the joints at least once every six months, or as often as the instrument maker decides in the patient's individual situation, for functionality and wear and tear.

Package content


- MultiMotion corrective joint
- Allen screwdriver
- MultiMotion Technical instructions
- MultiMotion Direction of force Quick Start Guide
- MultiMotion Patient documentation
- MultiMotion Tips & Tricks

The corrective joint is available in different variations. Check the package label to make sure the shown article number corresponds with the content of the package.


Components

1. Joint upper part
2. Central housing
3. Central housing screw
4. Power output indicator
5. Locking/unlocking lever
6. Joint lower part
7. Proximal modular bar adapter


8. Spring tension adjustment screw
9. Gear wheel
10. Distal modular bar adapter
11. Flexion / extension settings
12. Flexion / extension screws
13. Inclinement indicator
14. Bar attachment bushes


Dimensions


MM-100 / MM100I (Regular)


MM-200 / MM-200I (Small)


MM-110 / MM-110I (Regular)


MM-210 / MM-210I (Small)


MM-110US / MM-110IUS (Regular)

MM-210US / MM-210IUS (Small)


Product & Components (Regular)


No.	Art. No.	Description	Per
1A	MM-100	Corrective joint	Unit
1B	MM-100I	Corrective Joint with modular bar adapter	Unit
1C	MM-110	Corrective joint T-Bar	Unit
	MM-110US	Corrective joint T-Bar offset	Unit
1D	MM-110I	Corrective System Joint with integrated T-Bar	Unit
	MM-110IUS	Corrective System Joint, offset with integrated T-Bar	Unit
2	MM-1001	Central housing screw	Unit
3	MM-1002	Retainer ring	Unit
4	MM-1003	Housing	Unit
5	MM-1004	Force output indicator disc	Unit
6	MM-1005	Tension spring 10,2 Nm	Unit
7	MM-1006	Flexion / Extension stops	Set (2x)
8	MM-1007	Bar screws	Set (4x)
9	MM-1008	Bar attachment bushes	Set (4x)
10	MM-900	Allen screwdriver - 2.5mm	Unit

Product & Components (Small)

No.	Art. No.	Description	Per
1A	MM-200	Corrective joint	Unit
1B	MM-200I	Corrective Joint with modular bar adapter	Unit
1C	MM-210	Corrective joint T-Bar	Unit
	MM-210US	Corrective joint T-Bar offset	Unit
1D	MM-210I	Corrective System Joint with integrated T-Bar	Unit
	MM-210IUS	Corrective System Joint, offset with integrated T-Bar	Unit
2	MM-2001	Central housing screw	Unit
3	MM-2002	Retainer ring	Unit
4	MM-2003	Housing	Unit
5	MM-2004	Force output indicator disc	Unit
6	MM-2005	Tension spring 3,4 Nm	Unit
7	MM-2006	Flexion / Extension stops	Set (2x)
8	MM-2007	Bar screws	Set (4x)
9	MM-2008	Bar attachment bushes	Set (4x)
10	MM-900	Allen screwdriver - 2.5mm	Unit

Overview Bars & Accessories

No.	Art. No.	Description	Rigidity	Per
1	MM-701	Bar, straight 200 x 12 x 3.2 mm	Soft	Set (2x)
	MM-701B	Bar, straight 200 x 12 x 3.2 mm	Hard	Set (2x)
	MM-702	Bar, straight 200 x 16 x 4 mm	Soft	Set (2x)
	MM-702B	Bar, straight 200 x 16 x 4 mm	Hard	Set (2x)
2	MM-705	Bar, offset 350 x 12 x 3.2 mm	Soft	Set (2x)
	MM-705B	Bar, offset 350 x 12 x 3.2 mm	Hard	Set (2x)
	MM-706	Bar, offset 450 x 16 x 4 mm	Soft	Set (2x)
	MM-706B	Bar, offset 450 x 16 x 4 mm	Hard	Set (2x)
3	MM-709	Bar, T-Stirrup 150 x 12 x 3.2 mm	Soft	Set (2x)
	MM-709B	Bar, T-Stirrup 150 x 12 x 3.2 mm	Hard	Set (2x)
	MM-710	Bar, T-Stirrup 150 x 16 x 3.2 mm	Soft	Set (2x)
	MM-710B	Bar, T-Stirrup 150 x 16 x 3.2 mm	Hard	Set (2x)
4	MM-950	Ausrichtwerkzeug		Set

Joint selection

Determine which corrective joint unit is suitable to support the patient's joint. The MultiMotion "Small" model has a maximum torque output of 3.4 Nm. The "Regular" variation has a maximum torque output of 10.2 Nm.

MultiMotion joints are suitable for both children and adults. Both models are adjustable in angles of 0° to 120°. The joint is supplied ready for assembly including an allen screwdriver for setting the direction of force, the spring tension and to set the flexion and extension stops.

Note: make sure the joint is properly mounted, meaning the MultiMotion-logo on the housing should be readable (not upside down).

The joints can be used at both the left and right side and can be set with force in flexion or extension direction. This direction of force depends on the side of the body and the patient's joint for which the joint unit is used (see cover, Fig. A). The tension spring is painted red on one side and yellow on the other. The quick start guide and these instructions indicate which color at which side of the body shows the joint-specific direction of force. The standard assembly of the tension spring is with the red side facing upward.

Joint selection

Body Joint	Body weight / height	Corrective Joint		Free motion joint/R.O.M. joint	
		Small	Regular	Small	Regular
Wrist	< 22 kg	X		X	
	> 22 kg	X		X	
Elbow	< 22 kg	X		X	
	> 22 kg		X		X
Knee	< 22 kg	X		X	
	> 22 kg		X		X
Ankle	< 22 kg	X		X	
	> 22 kg		X		X
Hip (Children only)	< 105 cm	X			
	> 105 cm		X		

Selection flexion / extention

Body Joint	Body side	Flexion force		Extension force	
		Yellow	Red	Yellow	Red
Wrist	left	X			X
	right		X	X	
Elbow	left		X	X	
	right	X			X
Knee	left	X			X
	right		X	X	
Ankle	left	X			X
	right		X	X	

Changing the direction of force

The following description assumes that the joint is delivered ex-factory and therefore the tension spring has been assembled with the red side facing up. Make sure the joint faces with the logo readable in front of you, is in the neutral position (0°) as to minimize spring force, and the joint is unlocked (Fig. 1).

Step 1

Loosen the central screw on the housing using the allen screwdriver. Remove the screw and the underlying retainer ring (Fig. 2).

Step 2

Lift the housing over the central axis and lay it aside. Now the entire power indicator disc with the red indicator and circles is visible. Lift this also from the axis and lay it aside safely. Using your fingers, move the spring up over the central joint axis and remove it (Fig. 3).

Step 3

Now turn the spring around its axis with the yellow side facing upward, flex the joint 90° allowing you to replace the spring over the central axis. CAUTION: Make sure the spring is placed correctly between the recess of the central axis (Fig. 4).

Step 4

Using a cloth, wipe the power indicator disk clean and place it – but now with the yellow indicator and circles facing upward – back over the central axis. CAUTION: The tail of the spring should be positioned in the distal recess of the joint (Fig. 5).

Step 5

Place the housing and retainer ring back correctly and screw them back onto the joint (Fig 6).

Aligning the joint

In order for the joints to function correctly, it is important they are aligned correctly to one another, and a construction with sufficient torsion rigidity is built by using the free motion joint or R.O.M. joint. Make sure to take the necessary precautions, or use the alignment device (Art. Nr. MM-950).

Bending the joint

The material from both the corrective joint and the free motion / R.O.M. Joint is proximal and distal partly depleted (just below and above the Bar adapter). This allows to bend the joint as needed up to a maximum angle of 30°. Do not use bending irons but rather position the joint with a mounted Bar horizontally in a vice and push / pull the hinge into the desired position by hand.

Mounting the joint onto the Bars

Various Bar models are available in several widths and thicknesses. Select the Bar you need. Once you have adjusted the Bars to the contours of the limb, place the bars with the countersunk part facing down (!) over the attachment bushes of the joint. Use the supplied screws to fasten the bars and secure them with Loctite. TIP: If needed, you may loosen the attachment bushes and place them the other way around, allowing you to mount the Bars underneath instead of on top.

Locking the joint

In order to more easily position and remove the orthosis, the joint can be locked at virtually each desired angle. Distally from the power indicator window is the locking mechanism. The locking catch points distal as long as the joint is unlocked. In order to lock the joint, the catch should be pushed forward/up.

CAUTION

BEWARE OF RISK OF INJURY BY UNLOCKING. BEFORE UNLOCKING THE JOINT, NOTE THE SPRING TENSION SETTING. IF TENSION IS SET HIGH, IT MAY BE ADVISABLE TO HAVE A CARE GIVER ASSIST WITH UNLOCKING IN ORDER TO GRADUALLY GUIDE THE ORTHOSIS TOWARDS THE SET FLEXION OR EXTENSION FORCE.

Setting flexion and extension stops

There are two screws at the back of the corrective MultiMotion joint to set the range of motion. The standard stop setting in both corrective joints is 90° flexion and 90° extension. In the "Small" model, these flexion and extension stops can be set in 14° increments. In the "Regular" model, these increments are 12°. Use the supplied allen screwdriver to loosen the flexion and extension screws, position them as needed and retighten them.

Setting the spring tension

The allen screw at both sides of the joint is to set the spring tension using the supplied allen screwdriver. In each phase, the optimal objective is to realize a wearing time of 7-8 hours a day, uninterrupted without complications.

1st phase

Adaptation phase: set spring tension at minimal tension (Fig. A) and put the orthosis on the patient (Fig. B).

2nd phase

Work up to sub-maximum spring tension: this is the maximum spring tension the patient can tolerate for a period of 7-8 hours a day. Finding the sub-maximum spring tension is done progressively.

This means: start at the lowest tension and gradually increase the spring tension (Fig. C) until spasm is triggered. At that point, the tension should be slightly reduced (Fig. D) in order to relieve the spasm.


3rd phase

Rehabilitation: this is the phase where the patient's range of motion is progressively improved. As long as motion continues to be gained (Fig. E), spring tension is not increased. Only if motion is no longer gained, spring tension is increased (Fig. F) provided a wearing time of 7-8 hours a day

continues to be respected; this is normally done every 3-4 weeks. This phase is repeated until the patient's rehabilitation objective is achieved. With each change of force, prevention to trigger spasm should be a specific caution.

4th phase

Follow-up: this is the phase where the achieved objective (Fig. G) is retained. It is recommended to wear the orthosis daily for one hour in order to avoid recurrence of contractures (Fig. H).

Documenting Spring Tension

1st Step

The course of the LLPS therapy can best be observed by measuring of the Range of Motion. After sub-maximum spring tension is set (Fig. C and D), the next step is to determine the Active Range of Motion (A-ROM) while wearing the orthosis. This can be measured using a goniometer (Fig. I). The obtained result is to be documented in the appropriate chart (Fig. J).

2nd Step

Remove the orthosis from the patient. Documenting spring tension always has to be done with the joint locked and the orthosis in full extension (0° -position). If extension or flexion stops are used and prevent full extension just extend the orthosis up to the stop and lock the joint. The value shown on the indicator disc (Fig. K) is recorded in the chart (Fig. L).

At the next adjustment of the spring tension, the orthosis is again locked in full extension, and the set value recorded in the chart.

CAUTION





WHEN THE ORTHOSIS IS WORN BY THE PATIENT, THE INDICATOR DISC ROTATES ALONG WITH THE ORTHOSIS AND THUS, AT THAT MOMENT DOES NOT INDICATE ANY VALUE OF SPRING TENSION SETTINGS.

While the patient is wearing the orthosis the indicator disc in the joint turns along with the orthosis. Consequently, at that time it gives no indication of the set spring tension.

3rd Step

Measure the Passive Range of Motion (P-ROM) (Fig. M), again using a goniometer, without the patient wearing the orthosis and record the result in the appropriate chart (Fig. N).

Spring tension will be increased, as described in the 3rd stage - the rehabilitation phase - on several occasions. To keep track of the progress of the treatment, the results are measured and listed.

Cleaning instructions

The synthetic joint components may be cleaned with naphtha, but are unsuitable for aggressive agents such as acetone or similar. Try to keep the joint free from dust and dirt as much as possible. The joint should always be sufficiently greased. Use a standard white grease (e.g. Molycoat GS4500) to grease the joint.

Liability

Basko Healthcare as the manufacturer is solely liable if the product is used in accordance with the prescribed conditions and the purpose for which it is intended. Basko Healthcare recommends handling the product in accordance with applicable rules and performing maintenance as described in these user instructions. Non-compliance with the above may cause the MultiMotion joints to not function or function inadequately and may impact guarantee.

Changes, adjustments and modifications to the product not described in these instructions, are the responsibility and liability of the person who performs them.

CE conformity

These products meet requirements in Directive 93/42/EEC for medical devices. These products have been included in Class I based on classification criteria for medical devices as per annex IX of the Directive. Basko Healthcare therefore wrote up this conformity statement entirely at its own responsibility as per annex VII of the Directive.

Basko Healthcare

Pieter Lief tinckweg 16 | 1505 HX Zaandam
Tel.: +31 (0) 75 - 613 15 13 | Fax: +31 (0) 75 - 612 63 73
E-mail: verkoop@basko.com | Internet: www.basko.com

basko.com