# MultiMotion

Dynamic Hip abduction System

Patient documentation

**Basko**  
Healthcare  
basko.com

MM-9985 / Stand: 12/2013

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
## UK MultiMotion - Setting Spring Tension

The screw at the side of the joint is to set spring tension using the supplied roundhead screwdriver. In each phase, the optimal objective is to realize a wearing time of 7-8 hours a day, uninterrupted without counter response.

### 1<sup>st</sup> Phase

**Adaption phase:** set the spring tension at minimum force (Fig. 1) and apply the orthosis to the patient (Fig. 2).

### 2<sup>nd</sup> Phase

**Work up to sub-maximum spring tension:** this is the maximum spring tension the patient can tolerate for a time period of 7-8 hours a day. Finding the sub-maximum spring tension is done progressively. This means: start at the lowest tension and gradually increase spring tension (Fig. 3) until spasm is triggered. At that point, the tension should be slightly reduced in order to relieve the spasm (Fig. 4).

### 3<sup>rd</sup> Phase

**Rehabilitation:** this is the phase where the patient's range of motion is progressively improved. As long as motion continues to be gained, the spring tension is not increased (Fig. 5). Only if motion is no longer gained, the spring tension is increased (Fig. 6), provided a wearing time of 7-8 hours a day continues to be respected; this is normally done every 3-4 weeks. This phase is repeated until the patient's rehabilitation objective is achieved. With each change of force, prevention to trigger spasm should be avoided.

### 4<sup>th</sup> Phase

**Follow-up:** this is the phase where the achieved objective is retained (Fig. 7). It is recommended to wear the orthosis daily for one hour in order to avoid recurrence of contractures (Fig. 8).

## UK MultiMotion - Documenting Spring Tension

### 1<sup>st</sup> Step

The course of the LLPS therapy can best be observed by measuring of the Range of Motion. After sub-maximum spring tension is set (Fig. 3 and 4), the next step is to determine the Active Range of Motion (A-ROM) while **wearing** the orthosis. This can be measured using a goniometer (Fig. 9). The obtained result is to be documented in the appropriate chart (Fig. 10).

### 2<sup>nd</sup> Step

Remove the orthosis from the patient. Documenting spring tension always has to be done with the joint **locked** and the orthosis in **full extension** (0°-position). If extension or flexion stops are used and prevent full extension just extend the orthosis up to the stop and lock the joint. The value shown on the indicator disc (Fig. 11) is recorded in the chart (Fig. 12). At the next adjustment of the spring tension, the orthosis is again locked in full extension, and the set value recorded in the chart.

⚠ While the patient is wearing the orthosis the indicator disc in the joint turns along with the orthosis. Consequently, at that time it gives no indication of the set spring tension.

### 3<sup>rd</sup> Step

Measure the Passive Range of Motion (P-ROM) (Fig. 13), again using a goniometer, **without** the patient wearing the orthosis and record the result in the appropriate chart (Fig. 14).

Spring tension will be increased, as described in the 3rd stage - the rehabilitation phase - on several occasions. To keep track of the progress of the treatment, the results are measured and listed.