NL

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Material Composition: 90% Homopolymer Polypropylene, 10% Velcro hook and loop.

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Do not Iron

ENG

Purpose of Product: The Posterior Leaf Spring/Ankle Foot Orthosis provides assistive support and comfort while allowing patient to perform dynamic, functional activities. Fabricated in a predorsiflexed position to provide superior dorsiflexion assist.

Indications:

- Mild to moderate neurological deficit resulting from drop foot
- Flaccid foot-drop
- Peroneal Palsy

Contraindications:

- Broken skin or superficial wounds
- Allergy to component

Side Effects: No known side effects with correct use.

Fitting Instructions:

1. Apply a full-length cotton sock to patient's leg.
2. Insert brace in patient's standard shoe.
3. Insert foot into shoe. Wrap strap around calf and tighten firmly to create equal anterior-posterior compression.

Precautions: Follow the advice of the healthcare practitioner who prescribed or supplied the product. Wearing this product, as with any immobilizer product, could cause irritation (redness). Should this occur, readjust, or remove the product and contact your healthcare practitioner. Do not re-use the product for another patient. Any serious incident that occurs in connection with the product must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient are/is domiciled.

Care Instructions:

- Clean orthosis with mild soap and water. Towel dry.

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